



**Best Prices.  
Best Kits.  
Guaranteed.**

| <b>Drug</b> | <b>Cutoff Level</b> | <b>Approximate Detection Time in Urine</b> |
|-------------|---------------------|--|
| THC         | 50 ng/ml            | 0-30 Days                                  |
| COC         | 300 ng/ml           | 2-4 Days                                   |
| METH        | 300 ng/ml           | 3-5 Days                                   |
| OPI         | 300 ng/ml           | 2-4 Days                                   |
| AMP         | 300 ng/ml           | 2-4 Days                                   |
| BZO         | 200 ng/ml           | 3-7 Days                                   |

**NOTE:** The above chart gives approximate detection periods for each substance by test type. The ranges depend on amount and frequency of use, metabolic rate, body mass, age, overall health, drug tolerance, and urine pH.